

# BruxZir® NOW

“The Strength You Trust”



## The **only fully sintered zirconia** for in-office use

- Recommended for posterior and anterior crowns demanding maximum strength
- Mills in under 40 minutes\*
- Available exclusively for Glidewell milling systems
- The No. 1 brand of full-contour zirconia, with over 20 million restorations prescribed worldwide



MAXIMUM  
STRENGTH



NO OVEN  
NEEDED



MINIMAL  
PREPARATION

Before



After

This patient presented with failing restorations and recurrent caries on teeth #12-15. The replacement BruxZir® NOW restorations are esthetic, comfortable and strong enough to last the patient long into the future.

“

*I can seat BruxZir NOW crowns in 45 minutes without messing with ovens or glazing. Simple, easy and accurate, meaning I spend my time with my patients, not acting as a lab tech.*

”

– Maurice Trusas, DDS | Redding, CA

For more information  
[glidewelldental.com](http://glidewelldental.com) ■ 844-949-7184



## BENEFITS

- Select from regular and anterior formulations to meet the needs of each case
- Deliver restorations right from the mill with no post-processing required
- Easy delivery — conventional cementing with no need to bond
- Preserve your patients' existing dentition with minimal preparation requirements

## SPECIFICATIONS

### MATERIAL

- Pre-shaded, fully sintered zirconia
- Flexural strength: 1,100 MPa

### IDEAL FOR

- Single-unit anterior and posterior restorations
- Cases with limited interocclusal space
- Restorations with < 1 mm margin thickness

### SHADES

- 10 VITA® Classical Shades: A1, A2, A3, A3.5, B1, B2, C1, C2, D2, D3
- Posterior blocks also available in: A4, B3, B4, C3

## READY TO DELIVER IN UNDER 45 MINUTES



MILL

+



FINISH/POLISH

+




NO OVEN  
REQUIRED

=




TOTAL TIME

## EASY PREPARATION



- 1.25 mm ideal reduction (0.7 mm minimum)
- Chamfer margins preferred
- Labial and lingual walls must be convergent
- Labial surface should be cut in three planes



- 0.7 mm incisal/occlusal reduction required (1.5 mm recommended)
- Chamfer or shoulder margins preferred; feather-edge OK
- All axial walls must be convergent

\*Milling times vary based on prescription and mill used. Average milling time shown reflects use of the fastmill.io™ In-Office Mill.  
VITA is a registered trademark of VITA Zahnfabrik.

MKT-10317\_2 PK-268534-021820

For more information  
[glidewell.com](http://glidewell.com) ■ 844-949-7184

