


NEW!

Obsidian® **NOW***"Translucent Beauty, Timeless Strength"*

## The Lithium Disilicate with Easy Delivery and Great Results

- Recommended for anterior crowns
- Mills in under 35 minutes\*
- Available exclusively for Glidewell milling systems
- Excellent long-term clinical performance



PATIENT-PLEASING  
RESULTS



NO OVEN  
NEEDED



MINIMAL  
PREPARATION

### Obsidian® NOW Chairside Workflow



For more information  
[glidewelldental.com](http://glidewelldental.com) | 844-949-7184



# BENEFITS

- Deliver Obsidian® NOW right from the mill with no oven required
- Easy delivery — conventionally cemented or adhesively bonded
- Preserve your patients' existing dentition with minimal preparation requirements

# SPECIFICATIONS

## MATERIAL

- Versatile monolithic glass-ceramic
- Flexural strength: 476 MPa


## IDEAL FOR

- Single-unit anterior restorations
- Cases with limited interocclusal space
- Restorations with  $\geq 1.2$  mm thickness for circular and  $\geq 1.5$  mm for incisal/occlusal

## SHADES

- Classical Shades (A1, A2, A3, B1)


## READY TO DELIVER IN UNDER 45 MINUTES




# PREP GUIDE FOR ALL-CERAMIC RESTORATIONS

## ANTERIOR FULL-COVERAGE CROWN

### A. Full anterior wall thickness



### B. Shoulder



For more information  
[glidewelldental.com](http://glidewelldental.com) | 844-949-7184



\*Milling times vary based on prescription and mill used. Average milling time shown reflects use of the fastmill.io™ In-Office Mill. VITA is a registered trademark of VITA Zahnfabrik.

PK-5033462-060124