

# FASTPRINT.IO™: THE SIMPLE, INTUITIVE, HANDS-OFF 3D PRINTING SOLUTION

COMPETITOR

NEW!



fastprint.io™

3D PRINTING SOLUTION

Appliances from other printers are cloudy and use excessive material. The fastprint.io™ uses Digital Polishing™ for crystal-clear appliances with less material waste.


- **EXPANDED SAME-DAY SERVICES:**  
Print nightguards, retainers and other appliances in an easy, same-day workflow.
- **SIMPLE 1-CLICK DESIGNS AND DIGITAL POLISHING™:**  
Just one click, and you're ready to print crystal-clear appliances with minimal polishing needed.
- **AFFORDABILITY AND SAVINGS:**  
If you prescribe six or more nightguards per month, save over \$504 each month by printing them with the fastprint.io™ 3D Printing Solution instead of using a lab.\*

Plus, get direct support and training from Glidewell!




## Scan

Capture a digital impression of the patient's mouth using your preferred intraoral scanner.




## 1-click Design

Open fastdesign.io, fill out Rx, choose the appliance type, and let the AI create a design.



## Print

Press print on your fastprint.io to start printing.



## Post Process

Fully automated washing and curing stations allow for unattended operation. Just remove the supports, and you're done.



## Deliver





LEARN MORE  
[glidewell.io/fastprint](https://glidewell.io/fastprint)  
844-949-7184

\*Savings estimate based on comparison between \$90 lab price per nightguard versus material cost of 3D printing a nightguard.  
LuxCreo Digital Polishing is a trademark of LuxCreo Inc. Corporation.






 **Glidewell**

*for the sake of smiles*

# EXPAND YOUR SAME-VISIT CLINICAL SOLUTIONS

GLIDEWELL RESINS	TuffSplint™ Appliance Resin for Same-Day Nightguards	Quick and Affordable Modeling Resin
	 <p><b>\$3</b> resin cost per nightguard**</p>	 <p><b>\$3.50</b> resin cost per model**</p>

LUXCREO RESINS	DCA for Clear Aligners and Retainers	DNG for Clear Nightguards	DSG for Surgical Guides	DMR III for Dental Models
	 <p><b>\$6</b> resin cost per aligner or retainer*</p>	 <p><b>\$6</b> resin cost per nightguard*</p>	 <p><b>\$4</b> resin cost per guide**</p>	 <p><b>\$3</b> resin cost per model*</p>


3RD PARTY RESINS	detax Highend Medical Materials	RODIN® 3D Resin	Desktop Health (Flexcera™) Resins	keyprint® Resins	Dreve (Dentamid) Resins
	 <p>detax Freeprint® Crown: Permanent Crowns, Permanent Denture Teeth, and Long-Term Bridges</p>	 <p>Rodin Denture Base 2.0: Full and partial removable denture bases Rodin Titan: All-on-X and denture tooth arches Rodin Sculpture 2.0: Crowns, bridges, and inlays and onlays</p>	 <p>Flexcera Smile Ultra+: Permanent denture teeth, crowns, bridges, All-on-X, and inlays and onlays Flexcera Base Ultra+: Full and partial removable denture bases</p>	 <p>KeyOrtho IBT: Indirect Bonding Trays KeySplint Soft: Splints, Nightguards, and Bleaching Trays</p>	 <p>Dentamid FotoDent IBT: Indirect Bonding Trays Dentamid FotoDent tray2: Impression Trays (red, green and blue transparent) Dentamid FotoDent gingiva: Dental Gingiva Masks Dentamid FotoDent guide: Dental Surgical Guides</p>

Start printing with fastprint.io™ for only \$16,990\*

Includes washing and curing stations plus starter kit.

Contact us for our latest promotions and financing options!

\*Price does not include shipping or applicable taxes.



LEARN MORE  
[glidewell.io/fastprint](https://glidewell.io/fastprint)  
 844-949-7184

**Glidewell**  
*for the sake of smiles*

All third-party trademarks are property of their respective owners.

\*Cost estimates for Glidewell resins are based on internal data gathered using the fastprint.io with Glidewell TuffSplint Appliance Resin and Glidewell Modeling Resin. Cost estimates for LuxCuro resins are based on internal data gathered using the fastprint.io with DCA, DNG, DSG and DMR III.

PK-3887959-121024