

# fastscan.io™

## SCANNING SOLUTION

Includes a portable cart, laptop, and monitor.



Medit® i700®



Medit® i700® wireless



Medit® i900®



## Your Digital Workflow Starts Here

The fastscan.io™ Scanning Solution combines Medit® technology with fastscan.io software to deliver a complete, stress-free and affordable introduction to digital dentistry.

- Pay zero subscription or licensing fees.
- Get free **support and training from Glidewell.**
- Send digital cases to Glidewell for a **3-day turnaround and \$20 savings** plus reduced shipping costs, or send cases to your lab of choice.
- Choose your scanner bundle with the **Medit i700® wired, Medit i700 wireless, or the new Medit i900®.**

## What Your Peers Are Saying About Their Experience With fastscan.io



“

*fastscan.io lives up to the name — its speed and simplicity are exactly what we needed to expand our in-office treatments at a fraction of the cost of other scanners on the market. The fact that it has no monthly fees and even includes a cart, laptop, and monitor is outstanding.*”

– Joshua Prentice, DDS  
Flower Mound, TX



“

*Some in my staff had never used a scanner before. But after just two tries with fastscan.io, their scans looked just as good as mine! The training was great, the hardware is easy to use, and the Glidewell lab is wonderful. I can't say enough great things.*”

– George Marquis, DDS  
Fulton, MS



“

*I just can't believe the excellent customer service I get with the fastscan.io Scanning Solution — when I call with a question, they have an immediate answer. That's chairside support that I'm just so happy with. I feel like I have somebody in my corner.*”

– Nicole DeVecco Kim, DMD  
East Greenwich, RI



Scan the QR code for more information.

[glidewell.io/fastscan-io](https://glidewell.io/fastscan-io)

844-949-7184


 **Glidewell**

*for the sake of smiles*

## Three Solutions for Seamless Intraoral Scanning

	fastscan.io featuring Medit i700	fastscan.io Wireless featuring Medit i700	fastscan.io featuring Medit i900
<b>Imaging Technology</b>	3D-in-motion video technology/3D full color streaming capture	3D-in-motion video technology/3D full color streaming capture	3D-in-motion video technology/3D full color streaming capture
<b>Accuracy</b>	10.9 $\mu\text{m} \pm 0.98$	10.9 $\mu\text{m} \pm 0.98$	10.9 $\mu\text{m} \pm 0.98$
<b>Weight</b>	245 g	328 g (including battery, tip)	165 g
<b>Field of View</b>	15 x 13 mm	15 x 13 mm	L: 18 x 15 mm; M: 14 x 11.5 mm; S: 10 x 8 mm (sold separately)
<b>Adaptive Anti-Fogging</b>	✓	✓	✓
<b>Autoclavable Tip</b>	Up to 100 times	Up to 100 times	Up to 150 times
<b>UV-C Disinfection</b>	✓	✓	✓
<b>Remote Control Mode</b>	✓	✓	✓
<b>Free Training and Support</b>	✓	✓	✓
<b>Includes Cart, Laptop and Monitor</b>	✓	✓	✓
<b>Price</b>	\$14,990	\$18,500	\$20,999

**fastscan.io**<sup>TM</sup>  
SCANNING SOLUTION



**Maximize Your Scanning Potential With a fastscan.io<sup>TM</sup> Today!**

Contact us for our current promotional trade-in offers and payment plans on the fastscan.io<sup>TM</sup>.



If you prescribe as little as seven crowns per month, fastscan.io practically pays for itself. Scan the code to calculate your ROI and see how much you can save.  
[glidewell.io/fastscan-io-scanning-solution/roi-calculator](https://glidewell.io/fastscan-io-scanning-solution/roi-calculator)  
844-949-7184

 **Glidewell**  
*for the sake of smiles*