

NEW!

# BruxZir® NOW SRC

*“Implant Delivery Simplified”*



## The Fully Sintered Screw-Retained Zirconia Implant Crown for In-Office Use

- Recommended for implant cases demanding maximum strength
- Mills in under 40 minutes\*
- Available exclusively for Glidewell milling systems
- Facilitates delivery of same-day implant crowns
- Compatible with the Glidewell HT™ Implant System and NobelActive® implants from Nobel Biocare



MAXIMUM  
STRENGTH



NO OVEN  
NEEDED



PRE-CEMENTED  
ABUTMENT

For more information  
[glidewelldental.com](http://glidewelldental.com) | 844-949-7184



## BENEFITS

- Deliver implant crowns right from the mill in a single appointment
- Pre-cemented titanium abutment allows for easy delivery
- CAD/CAM fabrication ensures esthetic gingival margins and a natural emergence profile

## SPECIFICATIONS

### MATERIAL

- Pre-shaded, fully sintered zirconia
- Pre-installed titanium base
- Flexural strength: >1,000 MPa

### IDEAL FOR

- Prosthetic rehabilitation

### SHADES

- VITA Classical Shades (A1, A2, A3, A3.5, B1, C2, D2)

## READY TO DELIVER IN UNDER 45 MINUTES



MILL

+



FINISH/POLISH

+




NO OVEN  
REQUIRED

=



TOTAL TIME

## INSTRUCTIONS FOR USE



- (1) Fully sintered zirconia superstructure
- (2) Central screw access hole
- (3) Pre-cemented titanium abutment
- (4) Plastic mandrel
- (5) Abutment screw

IMPLANT SYSTEM	IMPLANT Ø	PLATFORM Ø
Hahn Tapered Implant System	Ø3.5 mm Ø4.3 mm Ø5.0 mm Ø7.0 mm	Ø3.5/4.3 mm Ø3.5/4.3 mm Ø5.0 mm Ø7.0 mm
Inclusive Tapered Implant System	Ø3.2 mm Ø3.7 mm Ø4.2 mm Ø4.7 mm Ø5.2 mm	Ø3.0 mm Ø3.5 mm Ø3.5 mm Ø4.5 mm Ø4.5 mm

BruxZir NOW SRC is compatible with the implant systems listed in the table above

For more information  
[glidewelldental.com](http://glidewelldental.com) | 844-949-7184



*for the sake of smiles*

\*Milling times vary based on prescription and mill used. Average milling time shown reflects use of the fastmill.io™ In-Office Mill.  
 VITA is a registered trademark of VITA Zahnfabrik.  
 The Glidewell HT Implant System is manufactured by Prismatic Dentalcraft, Inc., a wholly owned subsidiary of Glidewell Laboratories.  
 NobelActive is a registered trademark of Nobel Biocare Services AG. Glidewell HT is a trademark of Prismatic Dentalcraft, Inc.

PK-5158208-080124