



**With the glidewell.io<sup>TM</sup> In-Office Solution, effortlessly mill the crowns you want, when you want, with these special features:**

- Mill fully sintered BruxZir<sup>®</sup> crowns and bridges with no oven time required — just polish and cement, and you're done.
- Let MarginAI<sup>™</sup> and CrownAI<sup>™</sup> do the design work for you so you can spend your time with your patients.
- You're connected directly to Glidewell for training and support at any time.
- Love your scanner? Great! — glidewell.io is an open system and works with your preferred scanning system.

**“The glidewell.io In-Office Solution brings the ‘wow’ factor. In just one visit, your patient is receiving a beautiful BruxZir NOW crown that’s going to last a lifetime.”**

– Christina Do, DDS  
Costa Mesa, CA



**“With glidewell.io, the BruxZir NOW crowns are ready to deliver in 40 minutes or less. I don’t have to worry about oven-firing and then glazing, which allows me to send even more patients home in one quick visit.”**

– Lovedeep Randhawa, DDS  
Surrey, British Columbia



**“The glidewell.io system and its incorporated software are AMAZING! Scanning the patient and designing the restoration are easy to learn. The blocks mill quickly and come out of the mill basically polished and ready to go — CEREC<sup>®</sup> could never match that!”**

– Robert D. Rutner, DDS, DICOI  
Mountain View, CA



### FASTMILL.IO<sup>™</sup> PACKAGE

**Includes:**

- fastmill.io<sup>™</sup> In-Office Mill
- fastdesign.io<sup>™</sup> Software and Design Station
- Mill Accessory Kit
- One-Year Standard Warranty
- One-Year fastdesign.io Software License
- 45 Free Milling Blocks\*
- Free glidewell.io<sup>™</sup> Level I Advanced Course at Newport Beach, CA

**MSRP\*\* \$54,995**

*\*15 free milling blocks per month (any combination of blocks up to \$650 in value).*







*\*\*Prices do not include shipping or applicable taxes.*

*Opinions expressed in testimonials are solely the opinion of the doctor referenced and not those of Glidewell or its representatives.*

*CEREC is a registered trademark of Dentsply Sirona, Inc.*

# MILLING MATERIALS

The glidewell.io™ In-Office Solution offers dentists a versatile range of restorative materials designed to meet the specific needs of any clinical situation, whether for a full-coverage crown or a bridge. The fastmill.io™ In-Office Mill stands out as the only machine capable of producing same-day BruxZir® NOW—our fully sintered, oven-free restoration material.

	<b>BruxZir® NOW</b> 	<b>BruxZir® <i>esthetic</i> NOW</b> 	<b>CAMouflage® NOW</b> 	<b>BioTemps® NOW</b> 	<b>Obsidian® MILLING BLOCKS</b> 	<b>BruxZir® NOW SRC</b> 
<b>Materials</b>	Monolithic zirconia	Pre-shaded, fully sintered zirconia	Ceramic-filled nanohybrid resin	PMMA – Poly(methyl methacrylate)	Monolithic lithium silicate	Pre-shaded, fully sintered zirconia
<b>Indications</b>	Posterior crowns demanding maximum strength	Anterior crowns demanding high strength and maximum esthetics	Ideal for inlays, onlays and full-coverage crowns	Provisional crowns and bridges	Anterior and posterior crowns, veneers, inlays and onlays	Implant cases demanding maximum strength in the posterior
<b>Flexural Strength</b>	> 1,000 MPa	870 MPa	193 MPa	≥ 50 MPa	385 MPa	> 1,000 MPa
<b>Milling Time</b>	< 35 minutes for single-unit crowns < 2.5 hours for 3-unit bridges	< 35 minutes	< 10 minutes	< 10 minutes for single-unit crowns < 26 minutes for 3-unit bridges	< 20 minutes	< 40 minutes



**LEARN MORE**  
[glidewell.io/fastmillio-in-office-mill](http://glidewell.io/fastmillio-in-office-mill)  
 844-949-7184

**LEARN MORE ABOUT  
 OUR BRUXZIR CEMENT**   
[glidewelldirect.com/products/bruxzir-dual-cure-resin-cement](http://glidewelldirect.com/products/bruxzir-dual-cure-resin-cement)

