

BruxZir[®]
esthetic **NOW**

POSTERIOR BLOCK

“Beauty in a Single Visit”



NEW!



The Fully Sintered, Highly Esthetic Zirconia for In-Office Use

- Recommended for posterior crowns demanding high strength and excellent esthetics
- Mills in under 40 minutes*
- Available exclusively for Glidewell milling systems
- Optimized combination of strength and beauty



PATIENT-PLEASING
RESULTS



NO OVEN
NEEDED



POSTERIOR
CROWNS



HIGH
STRENGTH

For more information
glidewell.com/milling-blocks | 844-949-7184

glidewell.io[™]
IN-OFFICE SOLUTION

BENEFITS

- Lifelike translucency engineered for the posterior
- Deliver restorations right from the mill with no post-processing required
- Easy delivery — conventional cementing with no need to bond
- Preserve your patients' existing dentition with minimal preparation requirements

SPECIFICATIONS

MATERIAL

- Pre-shaded, fully sintered zirconia
- Flexural strength: 870 MPa


IDEAL FOR

- Single-unit posterior restorations, including bicuspid
- Cases with limited interocclusal space
- Restorations with < 1 mm margin thickness


SHADES

- 16 VITA® Classical Shades: A1, A2, A3

READY TO DELIVER IN UNDER 45 MINUTES



EASY PREPARATION



- 0.7 mm incisal/occlusal reduction required (1.5 mm recommended)
- Chamfer or shoulder margins preferred; feather-edge OK
- All axial walls must be convergent

For more information
glidewell.com/milling-blocks | 844-949-7184

*Milling times vary based on prescription and mill used. Average milling time shown reflects use of the fastmill.io™ In-Office Mill. VITA is a registered trademark of VITA Zahnfabrik.

 **Glidewell**
for the sake of smiles

PK-1524218-052025