

UNLOCK PRECISION WITH EASE YOUR WAY, WITH GLIDEWELL AND MEDIT

MEDIT *i*900 Mobility

MEDIT *i*900
CLASSIC



Scan the code to provide
your information and
we'll be in touch soon.


glidewell.io/fastscan-io
844-949-7184

 **Glidewell**
for the sake of smiles

YOUR DIGITAL WORKFLOW STARTS HERE

The fastscan.io™ Scanning Solution combines Medit® technology with Glidewell training and support to deliver a complete, stress-free and affordable introduction to digital dentistry.

- **Pay zero subscription or licensing fees.**
- Includes a portable cart, laptop, and monitor, **plus free support and training from Glidewell.**
- Send BruxZir® Zirconia model-free restorations cases to Glidewell for a **3-day turnaround and \$20 savings** or send cases to your lab of choice.
- Choose your scanner bundle with the **Medit i700 wired, Medit i700 wireless, Medit i900, Medit i900 Classic or Medit i900 Mobility.**



What Your Peers Are Saying About Their Experience With fastscan.io



"Glidewell has really nailed it—I like how it comes with the laptop, battery pack, a nice big monitor, and it's easy to move from room to room. **I would definitely recommend it to other practices looking to improve their digital workflow"**

— Kenneth Hsieh, DDS
Houston, TX



"Some in my staff had never used a scanner before. But after just two tries with fastscan.io, their scans looked just as good as mine! **The training was great, the hardware is easy to use, and the Glidewell lab is wonderful. I can't say enough great things."**



— George Marquis, DDS
Fulton, MS



"I just can't believe the excellent customer service I get with the fastscan.io Scanning Solution — when I call with a question, they have an immediate answer. That's chairside support that I'm just so happy with. I feel like I have somebody in my corner."

— Nicole DeVecco Kim, DMD
East Greenwich, RI

MEDIT i900 MOBILITY BEYOND WIRELESS



iPad device not included. Recommended iPad model: M4 (11 inch or 13 inch)

With Medit Link Express on iPad, Build Your Digital Dentistry Journey

1st iPad-Native Standalone App

A dedicated app built exclusively for iPad, optimized for mobile intraoral scanning workflows.

Optimized for Scan and Send

Focused workflow from patient registration to scanning and sending, enabling fast clinical turnaround.

Simple and Intuitive User Experience

Streamlined interface with only essential functions, enabling quick adoption with minimal learning curve.



Enjoy Direct Support and Training from Glidewell






Every fastscan.io purchase includes **free comprehensive training**, along with ongoing support. Our technicians are available weekdays from **5 a.m. to 5 p.m. Pacific** via phone, text, chat, or remote access.

INTERESTED IN TAKING YOUR DIGITAL JOURNEY FURTHER?

The fastscan.io™ Scanning Solution seamlessly integrates with glidewell.io™ In-Office Solution and fastprint.io™ 3D Printing Solution.

- Mill single-visit BruxZir crowns with no oven needed with fastmill.io™
- 3D print your own appliances with fastprint.io™
- Contact us for special offers with the fastscan.io



fastscan.io Scanning Solution bundle featuring Medit:					
Category	i900 mobility	i900 classic	i900	i700 wireless	i700
Model					
Scanning Technology	3D-in-motion video technology / 3D full color streaming capture				
Connection Type	Wireless	Wired		Wireless	Wired
Interface Type	Button		Touch	Button	
Dimensions	290 x 42 x 42 mm	223 x 37 x 35 mm		313 x 44 x 47 mm	248 x 44 x 47 mm
Total Weight	321 g	165 g		328 g	245 g
FOV (Scan Area)	18 x 15 mm			15 x 13 mm	
Camera	10 bit			8 bit	
Wireless Module	Wi-Fi 6E Module Built in, Bluetooth	N/A		60GHz Wireless USB 3.1Gen1 (hub) Bluetooth	N/A

Glidewell Exclusive:					
Free unlimited training and Support	Yes				
Includes cart, laptop and monitor	Yes				
Bundle Price	\$23,999	\$20,999	\$20,999	\$18,500	\$14,990
Scanner Only Price	\$21,999	\$18,999	\$18,999	\$15,500	\$13,500



If you prescribe as little as seven crowns per month, fastscan.io practically pays for itself. Scan the code to calculate your ROI and see how much you can save

glidewell.io/fastscan-io-scanning-solution/roi-calculator

844-949-7184

 **Glidewell**
for the sake of smiles

PK-5657295-080125