


NEW!

Obsidian® **NOW***"Translucent Beauty, Timeless Strength"*

The Lithium Disilicate with Easy Delivery and Great Results


- Recommended for anterior crowns
- Mills in under 35 minutes*
- Available exclusively for Glidewell milling systems
- Excellent long-term clinical performance



PATIENT-PLEASING
RESULTS



NO OVEN
NEEDED



MINIMAL
PREPARATION

Obsidian® NOW Chairside Workflow



For more information
glidewelldental.com | 844-949-7184



BENEFITS

- Deliver Obsidian® NOW right from the mill with no oven required
- Easy delivery — conventionally cemented or adhesively bonded
- Preserve your patients' existing dentition with minimal preparation requirements

SPECIFICATIONS

MATERIAL

- Versatile monolithic glass-ceramic
- Flexural strength: 476 MPa

IDEAL FOR

- Single-unit anterior restorations
- Cases with limited interocclusal space
- Restorations with ≥ 1.2 mm thickness for circular and ≥ 1.5 mm for incisal/occlusal

SHADES

- Classical Shades (A1, A2, A3, B1)


READY TO DELIVER IN UNDER 45 MINUTES




PREP GUIDE FOR ALL-CERAMIC RESTORATIONS

ANTERIOR FULL-COVERAGE CROWN

A. Full anterior wall thickness



B. Shoulder



For more information
glidewelldental.com | 844-949-7184



*Milling times vary based on prescription and mill used. Average milling time shown reflects use of the fastmill.io™ In-Office Mill. VITA is a registered trademark of VITA Zahnfabrik.

PK-5033462-060124