

fastsmile.ioTM
3D FACIAL SCANNER

BRING ESTHETICS INTO FOCUS



GET IN TOUCH

glidewell.io/fastsmile-3d

844-949-7184

 **Glidewell**

for the sake of smiles

SEE THE SMILE BEFORE YOU CREATE IT

The **fastsmile.io™ 3D Facial Scanner** is a 3D scanning device that captures high-resolution images of a patient's face for improved diagnostics, communication, and treatment planning.

- Capture full-face 3D scans to enhance smile designs, esthetic treatment planning, and restorative outcomes.
- Combine facial scans with intraoral data for esthetics that align with lip line, midline, and facial symmetry.

Modern dentistry isn't just about teeth—it's about how the smile fits the face. The **fastsmile.io™ 3D Facial Scanner** gives clinicians a powerful new way to visualize facial data and bring patients into the planning conversation with clarity and confidence.

Designed for speed and simplicity, **fastsmile.io** captures detailed facial scans in minutes, helping you align esthetics and function like never before.



“The fastsmile.io has been a great addition for esthetic planning. Having a full 3D facial scan makes it easier to visualize how the restoration will look in harmony with the patient's face—not just their smile.”

– Taylor Manalili, DDS | Irvine, CA

ELEVATE YOUR CASE PLANNING

Patients don't always understand dental terminology—but they understand what they see. Facial scans provide a clear, visual way to discuss smile goals, esthetics, and treatment possibilities, making consultations more engaging and easier to follow.

The result? Clearer conversations, stronger patient engagement, and more confident planning.


WORK WITH THE LAB

Integration with Glidewell's cloud-based platform allows the combined datasets to be submitted directly to the lab for collaborative case planning. Whether for veneers, crowns, or full-arch restorations, the scanner enhances the ability to preview final esthetics before initiating treatment. Dentists can simulate outcomes, gain patient approval, and confidently execute chairside workflows that meet both functional and cosmetic goals.

BUILT FOR TODAY'S ESTHETIC-FOCUSED DENTISTRY

Whether you're planning cosmetic cases, restorative work, or comprehensive smile designs, fastsmile.io™ helps you visualize outcomes in harmony with the face—supporting thoughtful, patient-centered planning.

The fastsmile.io™ fits naturally into your existing digital workflow, helping you move from scan to discussion without slowing down your day.



SEE THE DIFFERENCE

Bring facial awareness into your digital workflow with the **fastsmile.io™ 3D Facial Scanner**—and transform the way you plan, present, and communicate smiles.

- **Quick capture** – High-resolution 3D facial scans in minutes
- **Non-invasive** – Comfortable for patients, no contact required
- **Simple workflow** – Designed for effortless use by the dental team
- **All training included** – Includes training and ongoing support with Glidewell



Scanner Specifications

Model	FS-1-2025
Power and IO	Thunderbolt 4 USB-C and USB-A 3.0 for power and data
Operating Temperature	104° F maximum operating temperature
Dimensions (L x W x H)	240 x 44 x 115 mm
Weight	1.8 lbs
Field of View	200 x 110 mm
Accuracy	50 µm or better
Output Format	PLY, STL

PC Hardware Requirements

CPU	Intel i5
CPU Speed	2.7 GHz
GPU	NVIDIA GTX 1660 6 GB
Memory	32 GB RAM
Hard Drive Storage	500 GB NVMe
PC Connections	1x Thunderbolt 4 USB-C and 1x USB-A 3.0 port
Operating System	Windows 10 (64-bit) or later

This device is intended for visual planning and communication purposes only. It is not designed to diagnose, treat, cure, or prevent any disease or medical condition. Clinical decisions should be made solely by a licensed dental professional.



GET IN TOUCH

glidewell.io/fastsmile-3d
844-949-7184

